



More information on the website  
[mirror.radwag.com/tr/info,w1,QLH](http://mirror.radwag.com/tr/info,w1,QLH)

# 5 $\mu$ l ÷ 50 $\mu$ l - RW8-103-20-9 Otomatik Pipetler

PA-203-0003



The drawings, photos and graphics used are for illustrative purposes only.

## Teknik Veriler

### Metrolojik parametreler

Hacim	5 – 50 $\mu$ l
Systematic error	2,0% / 0,6%
Tekrarlanabilirlik error	2,0% / 0,3%
Uçlar	300 $\mu$ l



## İle uyumlu (Additional Fee)

XA 5Y.M.A.P pipetlerini kalibre etmek için Mikro Terazi  
MYA 5Y.P pipetlerini kalibre etmek için Mikro Terazi

Pipet kalibrasyonu için otomatik teraziler  
Pipet uçları