



More information on the website
mirror.radwag.com/it/info,w1,QPS

Cilindrico, set (1 g - 5 kg), imballaggi in plastica

OK-502-0010



The drawings, photos and graphics used are for illustrative purposes only.

Specifiche

| Parametri metrologici | |
|-----------------------|-----------------|
| Classe OIML | E1 |
| Costruzione | |
| Materiale | Stainless steel |
| Peso | 11 110 g |
| Valore nominale | 1 g – 5 kg |
| Numero | 16 pcs |

Compatibile con (Additional Fee)

Pinzette

Pennelli