



More information on the website  
[mirror.radwag.com/it/info,w1,QPS](https://mirror.radwag.com/it/info,w1,QPS)

# Cilindrico, set (1 g - 5 kg), imballaggi in plastica

OK-502-0010



The drawings, photos and graphics used are for illustrative purposes only.

## Specifiche

Parametri metrologici	
Classe OIML	E1
Costruzione	
Materiale	Stainless steel
Peso	11 110 g
Valore nominale	1 g – 5 kg
Numero	16 pcs

## Compatibile con (Additional Fee)

Pinzette

Pennelli