



More information on the website
mirror.radwag.com/en/info,w1,R3Z

Tips OMNITIP™ - Starter ST 10µl

PT-114-0001



The drawings, photos and graphics used are for illustrative purposes only.

Overview

OMNITIP™ Starter ST 10 µl tips are high quality laboratory accessories designed for precision pipetting. Designed with efficiency and convenience in mind, they are an indispensable tool for everyday laboratory work.

Capacity: 10 µl Colour: transparent Set contains: box, 6 racks and 576 tips

Application: - precise dispensing of biological samples, - accurate pipetting of medicines and medical solutions, - preparing and transferring chemical solutions, - an indispensable tool in experiments requiring precise dosing.

Compatible with (Additional Fee)

Workstation for Pipettes Calibration

Automatic Variable-Volume Pipettes