



More information on the website  
[mirror.radwag.com/en/info,w1,RDV](https://mirror.radwag.com/en/info,w1,RDV)

# M1 Mass Standard, Knob Weights, Set (1 g - 200 g), Plastic Box

OK-402-0066



The drawings, photos and graphics used are for illustrative purposes only.

## Datasheet

### Metrological parameters

OIML Class	M1
Weight	610 g
Nominal value	1 – 200 g
Quantity	11 pcs

## Compatible with (Additional Fee)

Brushes