



More information on the website
mirror.radwag.com/en/info,w1,RG1

F1 Mass Standard - Knob Weights Without Adjustment Chamber, Set (1 mg - 50 g), Plastic Box

OK-202-0082



The drawings, photos and graphics used are for illustrative purposes only.

Datasheet

Metrological parameters	
OIML Class	F1
Construction	
Material	Stainless steel
Weight	111.11 g
Nominal value	1 mg – 50 g
Quantity	20 pcs

Compatible with (Additional Fee)

WPY Multifunctional Scale
 TMX H Stainless steel Multifunctional Scale
 WLY Scales for PGC
 HY10 Multifunctional Scale
 PUE 5 Weighing Terminal
 HY10 H Stainless Steel Multifunctional Scale

WPY Scales for PGC
 HY10.H Scales for PGC
 Tweezers
 WLY Precision Balance
 Brushes

