



More information on the website  
[mirror.radwag.com/en/info,w1,RG1](https://mirror.radwag.com/en/info,w1,RG1)

# F1 Mass Standard - Knob Weights Without Adjustment Chamber, Set (1 mg - 50 g), Plastic Box

OK-202-0082



The drawings, photos and graphics used are for illustrative purposes only.

## Datasheet

Metrological parameters	
OIML Class	F1
Construction	
Material	Stainless steel
Weight	111.11 g
Nominal value	1 mg – 50 g
Quantity	20 pcs

## Compatible with (Additional Fee)

WPY Multifunctional Scale  
 TMX H Stainless steel Multifunctional Scale  
 WLY Scales for PGC  
 HY10 Multifunctional Scale  
 PUE 5 Weighing Terminal  
 HY10 H Stainless Steel Multifunctional Scale

WPY Scales for PGC  
 HY10.H Scales for PGC  
 Tweezers  
 WLY Precision Balance  
 Brushes

