



More information on the website
mirror.radwag.com/en/info,w1,RJ9

F1 Mass Standard - Knob Weights Without Adjustment Chamber, Set (1 g - 100 g), Wooden Box

OK-201-0062



The drawings, photos and graphics used are for illustrative purposes only.

Datasheet

| Metrological parameters | |
|-------------------------|-----------------|
| OIML Class | F1 |
| Construction | |
| Material | Stainless steel |
| Weight | 210 g |
| Nominal value | 1 – 100 g |
| Quantity | 9 pcs |

Compatible with (Additional Fee)

Tweezers

Brushes