



More information on the website  
[mirror.radwag.com/it/info,w1,RMM](http://mirror.radwag.com/it/info,w1,RMM)

# Pesi standard F1 - cilindrico - 1000 kg - cassa di legno

OP-201-0034



The drawings, photos and graphics used are for illustrative purposes only.

## Specifiche

| Parametri metrologici |                 |
|-----------------------|-----------------|
| Classe OIML           | F1              |
| Costruzione           |                 |
| Materiale             | Stainless steel |
| Valore nominale       | 1000 kg         |

## Compatibile con (Additional Fee)

Pennelli