



More information on the website  
[mirror.radwag.com/de/info,w1,RW7](https://mirror.radwag.com/de/info,w1,RW7)

# Knopfgewichte, Set (1 g - 2 kg), Kunststoff-Etui



The drawings, photos and graphics used are for illustrative purposes only.

## Technische Daten

| Messtechnische Parameter |            |
|--------------------------|------------|
| OIML-Klasse              | M1         |
| Gewicht                  | 6 110 g    |
| Nominalwert              | 1 g – 2 kg |
| Anzahl                   | 15 St.     |

## Kompatibel mit

Pinsel