



More information on the website  
[mirror.radwag.com/de/info,w1,RW7](http://mirror.radwag.com/de/info,w1,RW7)

# Prüfgewichte der M1 Klasse - Knopfgewichte, Set (1 g - 2 kg), Kunststoff-Etui

OK-402-0069



The drawings, photos and graphics used are for illustrative purposes only.

## Technische Daten

Messtechnische Parameter	
OIML-Klasse	M1
Gewicht	6 110 g
Nominalwert	1 g – 2 kg
Anzahl	15 St.

## Kompatibel mit (Additional Fee)

Pinsel