



More information on the website  
[mirror.radwag.com/en/info,w1,RW7](https://mirror.radwag.com/en/info,w1,RW7)

# M1 Mass Standard, Knob Weights, Set (1 g - 2 kg), Plastic Box

OK-402-0069



The drawings, photos and graphics used are for illustrative purposes only.

## Datasheet

### Metrological parameters

OIML Class	M1
Weight	6 110 g
Nominal value	1 g – 2 kg
Quantity	15 pcs

## Compatible with (Additional Fee)

Brushes