



More information on the website  
[mirror.radwag.com/it/info,w1,S4G](https://mirror.radwag.com/it/info,w1,S4G)

# Cilindrico, set (1 mg - 200 g), cassa di legno

OK-301-0039



The drawings, photos and graphics used are for illustrative purposes only.

## Specifiche

| Parametri metrologici |                 |
|-----------------------|-----------------|
| Classe OIML           | F2              |
| Costruzione           |                 |
| Materiale             | Stainless steel |
| Peso                  | 611,11 g        |
| Valore nominale       | 1 mg – 200 g    |
| Numero                | 23 pcs          |

## Compatibile con (Additional Fee)

Pinzette

Pennelli