



More information on the website
mirror.radwag.com/en/info,w1,SCH

Cable P0108

WX-013-0003



The drawings, photos and graphics used are for illustrative purposes only.

Overview

RS 232 cable: Balance - computer. Standard cable length: 1.7 m. Compatible with:

- Scales with PUE C315 terminal
- PUE C315 terminal

Compatible with (Additional Fee)

MA X7.IC.A Moisture Analyzer
WPY Multifunctional Scale
C32 Multifunctional Scale
MA 5Y Moisture Analyzer
C315 Platform scale
WLY Scales for PGC
AS X7 Analytical Balance
PS X7 Balances for PGC
C315.D Baby Scale
WLC Precision Balance
PS R2 Precision Balance
C315.4P2 Beam Scales
WLC X7 Precision Balance

C315.OW Personal Scale
XA 5Y.F Analytical Balance for Filters
WLY Precision Balance
UMA Automatic Mass Comparator
PS X7 Precision Balance
APP 5Y.KO Manual Mass Comparator
Robotic Mass Comparators
Density Measurement
WLC C/2 Precision Balance
MAS.R2 Modules
PS R1 Precision Balance
5Y.PM.KB Manual Mass Comparator
AS R2 PLUS Analytical Balance

C315.4 Platform Scale
C32.PM Precision Balance
MA X7.A Moisture Analyzer
C315.O Personal Scale
WPY Scales for PGC
C315.4P pallet scales
WAY 5Y.KO Manual Mass Comparator
MA 5Y.IC Moisture Analyzer
MPS.M.R2 Modules
MA R Moisture Analyzer
AS 5Y Analytical Balance
AK-4 Automatic Mass Comparator
PS 5Y Precision Balance
PS 5Y Balances for PGC
AKM-2 Automatic Mass Comparator
PS R2.H Precision Balance
WTC Precision Balance
H315 4C Platform Scale

PUE C32 Weighing Terminal
PUE 7.1 Weighing Terminal
Postal Scales for Packages
Automatic Vacuum Mass Comparators
PMV 5Y Moisture Analyzer
XA 5Y.KO Manual Mass Comparator
PUE C315 Weighing Terminal
MAS.R1 Modules
MPS.M.R1 Modules
UYA 5Y.KO Manual Mass Comparator
AS R1 PLUS Analytical Balance
C315.4B Bed Scale
MPS.T.R1 Modules
Nano-Comparator NANO.AK-4.500.5Y
C315.K Chair Scale
MAS.1.Y Modules
C315.8B Ramps Bed Scale