



More information on the website  
[mirror.radwag.com/de/info,w1,SKL](http://mirror.radwag.com/de/info,w1,SKL)

# Prüfgewichte der M1 Klasse - Knopfgewichte, Set (1 g - 50 g), Holz-Etui

OK-401-0071



The drawings, photos and graphics used are for illustrative purposes only.

## Technische Daten

### Messtechnische Parameter

OIML-Klasse	M1
Gewicht	110 g
Nominalwert	1 – 50 g
Anzahl	8 St.

## Kompatibel mit (Additional Fee)

Pinsel