



More information on the website  
[mirror.radwag.com/de/info,w1,TPQ](http://mirror.radwag.com/de/info,w1,TPQ)

# Prüfgewichte der E2 Klasse - Zylindergewichte, Satz (1 kg - 5 kg), Holz- Etui



The drawings, photos and graphics used are for illustrative purposes only.

## Technische Daten

Messtechnische Parameter	
OIML-Klasse	E2
Konstruktion	
Material	Edelstahl
Gewicht	10 000 g
Nominalwert	1 – 5 kg
Anzahl	4 St.