



More information on the website  
[mirror.radwag.com/de/info,w1,TUK](https://mirror.radwag.com/de/info,w1,TUK)

# Prüfgewicht der E1 Klasse - Plättchenform - 500 mg

OP-503-0045



The drawings, photos and graphics used are for illustrative purposes only.

## Technische Daten

### Messtechnische Parameter

OIML-Klasse E1

### Konstruktion

Material Edelstahl

Nominalwert 500 mg

## Kompatibel mit (Additional Fee)

Pinzetten

Pinself