



More information on the website
mirror.radwag.com/en/info,w1,U1A

SAL/H Stainless Steel Laboratory Antivibration Table

WX-001-0077



The drawings, photos and graphics used are for illustrative purposes only.

Datasheet

Construction	
Construction	Stainless steel AISI304
Load-bearing capacity	40 kg
Physical parameters	
Packaging dimensions W x D x H	820×711×816 mm
Net weight	61 kg
Gross weight	65-70 kg
Length	770 mm
Width	680 mm
Height	815±10 mm

Compatible with (Additional Fee)

AS X7 Analytical Balance
 PS R2 Precision Balance
 WLC X7 Precision Balance

WTC Precision Balance
 XA 5Y.F Analytical Balance for Filters
 WLC C/2 Precision Balance

WLC Precision Balance
C32.PM Precision Balance
PS X7 Precision Balance
PS R2.H Precision Balance

PS R1 Precision Balance
AS R2 PLUS Analytical Balance
AS R1 PLUS Analytical Balance
MAS.1.Y Modules

Device dimensions W x D x H

