



More information on the website  
[mirror.radwag.com/en/info,w1,U1A](https://mirror.radwag.com/en/info,w1,U1A)

# SAL/H Stainless Steel Laboratory Antivibration Table



The drawings, photos and graphics used are for illustrative purposes only.

## Datasheet

Construction	
Construction	Stainless steel AISI304
Load-bearing capacity	40 kg
Physical parameters	
Packaging dimensions	820×711×816 mm
Net weight	61 kg
Gross weight	65-70 kg
Length	770 mm
Width	680 mm
Height	815±10 mm

## Compatible with

MUYA Modules  
 MAS1 Modules  
 MAS Modules

MPS Modules  
 C32.PM Precision Balance  
 XA 5Y.F Analytical Balance for Filters

# Device dimensions

SAL/H

