



More information on the website  
[mirror.radwag.com/en/info,w1,U4E](https://mirror.radwag.com/en/info,w1,U4E)

## Suction pump for liquids

WX-001-0445



The drawings, photos and graphics used are for illustrative purposes only.

### Overview

Precision suction, error-free measurements The product is supplied with two internal pumps and two built-in buttons used to operate separate pumps. Backlit buttons inform you about operation of a specific pump. The set has two tips for sucking water out of a single beaker (one-channel tip) or four pipettes at a time (four-channel tip). Water is carried from both channels to the beaker through one tube. The power is provided to the pump through a USB A in the weighing terminal. You can supply power from the phone charger, for instance.

## Compatible with (Additional Fee)

Automatic Balances for Pipettes Calibration

XA 5Y.M.A Microbalance

XA 5Y.A Analytical Balance

XA 5Y Analytical Balance

XA 5Y.M Microbalance

MYA 5Y Microbalance