



More information on the website
mirror.radwag.com/it/info,w1,U53

Cilindrico, set (1 mg - 5 g), imballaggi in plastica

OK-102-0040



The drawings, photos and graphics used are for illustrative purposes only.

Specifiche

| Parametri metrologici | |
|-----------------------|-----------------|
| Classe OIML | E2 |
| Costruzione | |
| Materiale | Stainless steel |
| Peso | 6,11 g |
| Valore nominale | 1 mg – 5 g |
| Numero | 16 pcs |

Compatibile con (Additional Fee)

Pinzette

Pennelli