



More information on the website  
[mirror.radwag.com/it/info,w1,UAW](http://mirror.radwag.com/it/info,w1,UAW)

# Pipetta automatica 0,5 µl ÷ 10 µl - RW8-101-20-9

PA-203-0001



The drawings, photos and graphics used are for illustrative purposes only.

## Specifiche

Parametri metrologici	
Volume	0,5 – 10 µl
Systematic error	2,5% / 1,0%
Random error	1,5% / 0,8%
Tips	10 µl



## Compatibile con (Additional Fee)

Microbilancia per la calibrazione delle pipette XA 5Y.M.A.P  
 Microbilancia per la calibrazione delle pipette MYA 5Y.P

Puntali per pipette  
 Bilance automatiche per la calibrazione delle pipette