



More information on the website  
[mirror.radwag.com/de/info,w1,USF](http://mirror.radwag.com/de/info,w1,USF)

# Prüfgewichte der M1 Klasse - Knopfgewichte, Set (1 mg - 5 kg), Holz- Etui

OK-401-0085



The drawings, photos and graphics used are for illustrative purposes only.

## Technische Daten

Messtechnische Parameter	
OIML-Klasse	M1
Gewicht	11 111,11 g
Nominalwert	1 mg – 5 kg
Anzahl	28 St.

## Kompatibel mit (Additional Fee)

Pinsel