



More information on the website  
[mirror.radwag.com/it/info,w1,UYL](https://mirror.radwag.com/it/info,w1,UYL)

# Pesi standard F1 - cilindrico - 2000 kg - cassa di legno

OP-201-0035



The drawings, photos and graphics used are for illustrative purposes only.

## Specifiche

### Parametri metrologici

Classe OIML	F1
-------------	----

### Costruzione

Materiale	Stainless steel
-----------	-----------------

Valore nominale	2000 kg
-----------------	---------

## Compatibile con (Additional Fee)

Pennelli