



More information on the website  
[mirror.radwag.com/en/info,w1,VAG](https://mirror.radwag.com/en/info,w1,VAG)

# M1 Mass Standard, Knob Weights, Set (1 g - 500 g), Wooden Box

OK-401-0074



The drawings, photos and graphics used are for illustrative purposes only.

## Datasheet

Metrological parameters	
OIML Class	M1
Weight	1 110 g
Nominal value	1 – 500 g
Quantity	12 pcs

## Compatible with (Additional Fee)

Brushes