



More information on the website  
[mirror.radwag.com/it/info,w1,VAG](https://mirror.radwag.com/it/info,w1,VAG)

# Pesi standard M1, cilindrici, set (1 g - 500 g), cassa di legno

OK-401-0074



The drawings, photos and graphics used are for illustrative purposes only.

## Specifiche

### Parametri metrologici

Classe OIML	M1
Peso	1 110 g
Valore nominale	1 – 500 g
Numero	12 pcs

## Compatibile con (Additional Fee)

Pennelli