



More information on the website  
[mirror.radwag.com/en/info,w1,VBY](https://mirror.radwag.com/en/info,w1,VBY)

# F1 Mass Standard - Knob Weights Without Adjustment Chamber, Set (1 g - 200 g), Plastic Box

OK-202-0074



The drawings, photos and graphics used are for illustrative purposes only.

## Datasheet

Metrological parameters	
OIML Class	F1
Construction	
Material	Stainless steel
Weight	610 g
Nominal value	1 – 200 g
Quantity	11 pcs

## Compatible with (Additional Fee)

Tweezers

Brushes