



More information on the website  
[mirror.radwag.com/en/info,w1,VCQ](http://mirror.radwag.com/en/info,w1,VCQ)

# Knob weights, Set (1 g - 2 kg), Wooden Box

OK-501-0032



The drawings, photos and graphics used are for illustrative purposes only.

## Datasheet

Metrological parameters	
OIML Class	E1
Construction	
Material	Stainless steel
Weight	6 110 g
Nominal value	1 g – 2 kg
Quantity	15 pcs

## Compatible with (Additional Fee)

Tweezers

Brushes