



More information on the website  
[mirror.radwag.com/en/info,w1,VHG](https://mirror.radwag.com/en/info,w1,VHG)

# Knob weights, Set (1 kg - 5 kg), Plastic Box

OK-502-0012



The drawings, photos and graphics used are for illustrative purposes only.

## Datasheet

Metrological parameters	
OIML Class	E1
Construction	
Material	Stainless steel
Weight	10 000 g
Nominal value	1 – 5 kg
Quantity	4 pcs

## Compatible with (Additional Fee)

Brushes