



More information on the website  
[mirror.radwag.com/de/info,w1,VHN](http://mirror.radwag.com/de/info,w1,VHN)

# Prüfgewichte der M1 Klasse - Knopfgewichte, Set (1 mg - 50 g), Holz- Etui

OK-401-0079



The drawings, photos and graphics used are for illustrative purposes only.

## Technische Daten

Messtechnische Parameter	
OIML-Klasse	M1
Gewicht	111,11 g
Nominalwert	1 mg – 50 g
Anzahl	20 St.

## Kompatibel mit (Additional Fee)

Pinsel