



More information on the website  
[mirror.radwag.com/en/info,w1,VHN](https://mirror.radwag.com/en/info,w1,VHN)

# M1 Mass Standard, Knob Weights, Set (1 mg - 50 g), Wooden Box

OK-401-0079



The drawings, photos and graphics used are for illustrative purposes only.

## Datasheet

### Metrological parameters

|               |             |
|---------------|-------------|
| OIML Class    | M1          |
| Weight        | 111.11 g    |
| Nominal value | 1 mg – 50 g |
| Quantity      | 20 pcs      |

## Compatible with (Additional Fee)

Brushes