



More information on the website  
[mirror.radwag.com/en/info,w1,VXA](https://mirror.radwag.com/en/info,w1,VXA)

# Knob weights, Set (1 mg - 200 g), Plastic Box

OK-102-0028



The drawings, photos and graphics used are for illustrative purposes only.

## Datasheet

Metrological parameters	
OIML Class	E2
Construction	
Material	Stainless steel
Weight	611.11 g
Nominal value	1 mg – 200 g
Quantity	23 pcs

## Compatible with (Additional Fee)

Tweezers

Brushes