



More information on the website
mirror.radwag.com/de/info,w1,W1J

Prüfgewicht der E2 Klasse - Knopfgewicht - 500 g

OP-103-0047



The drawings, photos and graphics used are for illustrative purposes only.

Technische Daten

| Messtechnische Parameter | |
|--------------------------|-----------|
| OIML-Klasse | E2 |
| Konstruktion | |
| Material | Edelstahl |
| Nominalwert | 500 g |

Kompatibel mit (Additional Fee)

Pinsel

Pinzetten