



More information on the website
mirror.radwag.com/en/info,w1,W20

Under-pan Weighing Rack for AS and PS Balances



The drawings, photos and graphics used are for illustrative purposes only.

Overview

A rack for under-pan weighing intended for AS and PS balances.

The rack is designed for weighing loads requiring carrying out the measuring process outside the main weighing platform, e.g. loads with non-standard dimensions and sizes and those generating a magnetic field.

The rack also facilitates determining density of loads.

The rack is available in powder coated mild steel finishing.

Compatible with

AS X7 Analytical Balance
 PS X2 Precision Balance
 AS X2 PLUS Analytical Balances
 PS R2 Precision Balance
 PS X2 Balances for PGC
 WLC Precision Balance
 PS X7 Precision Balance

AS 5Y Analytical Balance
 PS R2.H Precision Balance
 WLC C/2 Precision Balance
 PS R1 Precision Balance
 AS R2 PLUS Analytical Balance
 PS 5Y Precision Balance
 AS R1 PLUS Analytical Balances

Device dimensions

