




More information on the website  
[mirror.radwag.com/us/info,w1,W20](https://mirror.radwag.com/us/info,w1,W20)

# Under-hook weighing rack for AS and PS balances

WX-002-0111



The drawings, photos and graphics used are for illustrative purposes only.

## Overview

A rack for under-hook weighing intended for balances AS and PS.

The rack is designed for weighing loads requiring carrying out the measuring process outside the main weighing platform, e.g. loads with non-standard dimensions and sizes and those generating a magnetic field.

The rack also facilitates determining density of loads.

The rack is available in powder coated mild steel finishing.

## Compatible with (Additional Fee)

AS X7 Analytical Balance  
PS X7 Balances for PGC  
PS R2 Precision Balance  
WLC Precision Balance  
AS 5Y Analytical Balance  
PS 5Y Precision Balance  
PS 5Y Balances for PGC

PS R2.H Precision Balance  
PS X7 Precision Balance  
WLC C/2 Precision balance  
PS R1 Precision Balance  
AS R2 PLUS Analytical Balance  
AS R1 PLUS Analytical Balance

## Device dimensions W x D x H

