

Under-hook weighing rack for AS and PS balances

WX-002-0111





The drawings, photos and graphics used are for illustrative purposes only.

Overview

A rack for under-hook weighing intended for balances AS and PS.

The rack is designed for weighing loads requiring carrying out the measuring process outside the main weighing platform, e.g. loads with non-standard dimensions and sizes and those generating a magnetic field. The rack also facilitates determining density of loads.

The rack is available in powder coated mild steel finishing.

Compatible with (Additional Fee)

AS X7 Analytical Balance PS X7 Balances for PGC PS R2 Precision Balance WLC Precision Balance PS X7 Precision Balance AS 5Y Analytical Balance PS 5Y Balances for PGC PS R2.H Precision Balance WLC C/2 Precision balance PS R1 Precision Balance AS R2 PLUS Analytical Balance PS 5Y Precision Balance AS R1 PLUS Analytical Balance

Device dimensions W x D x H





