



More information on the website
mirror.radwag.com/en/info,w1,W3Y

PSW/DUAL Professional Dual Weighing Workstation

WX-001-0402



The drawings, photos and graphics used are for illustrative purposes only.

Datasheet

Physical parameters	
Delivery components	table for calibration, anti-vibration table
Packaging dimensions	1140×820×986 mm
Net weight	220 kg (table 140 + stone 80)
Gross weight	252 kg
Length	1350 mm
Width	750 mm
Height	1100 mm
Number of workstations	2



Compatible with (Additional Fee)





MAS1 Modules
 XA 5Y.M.A.P Microbalance for pipettes calibration

AS 5Y Analytical Balance
 XA 5Y.F Analytical Balance for Filters

MYA 5Y.P Microbalance for pipettes calibration
 MYA 5Y.F Microbalance for Filters
 AS X7 Analytical Balance
 AS X2 PLUS Analytical Balances
 XA 5Y.M.A Microbalance
 MAS Modules
 MPS Modules

MYA 5Y Microbalance
 AS R2 PLUS Analytical Balance
 XA 5Y.M Microbalance
 UYA 5Y.KO Manual Mass Comparator
 AS R1 PLUS Analytical Balance
 XA 5Y.M.A.S Microbalance for stents weighing

Device dimensions



—SDKP

—SDKP DUAL