



More information on the website  
[mirror.radwag.com/en/info,w1,W3Y](http://mirror.radwag.com/en/info,w1,W3Y)

# PSW/DUAL Professional Dual Weighing Workstation

WX-001-0402



The drawings, photos and graphics used are for illustrative purposes only.

## Datasheet

Physical parameters	
Delivery components	table for calibration, anti-vibration table
Packaging dimensions W x D x H	1140×820×986 mm
Net weight	220 kg (table 140 + stone 80)
Gross weight	252 kg
Length	1350 mm
Width	750 mm
Height	1100 mm
Number of workstations	2



## Compatible with (Additional Fee)





XA 5Y.M.A.P Microbalance for pipettes calibration  
 MYA 5Y.P Microbalance for pipettes calibration

MYA 5Y Microbalance  
 AS R2 PLUS Analytical Balance

MYA 5Y.F Microbalance for Filters  
AS X7 Analytical Balance  
XA 5Y.M.A Microbalance  
AS 5Y Analytical Balance  
XA 5Y.F Analytical Balance for Filters

XA 5Y.M Microbalance  
UYA 5Y.K0 Manual Mass Comparator  
AS R1 PLUS Analytical Balance  
XA 5Y.M.A.S Microbalance for stents weighing  
MAS.1.Y Modules

## Device dimensions W x D x H



—SDKP

—SDKP DUAL