



More information on the website  
[mirror.radwag.com/en/info,w1,WRY](http://mirror.radwag.com/en/info,w1,WRY)

# F1 Mass Standard - Knob Weights With Adjustment Chamber, Set (1 g - 5 kg), Plastic Box

OK-202-0060



The drawings, photos and graphics used are for illustrative purposes only.

## Datasheet

Metrological parameters	
OIML Class	F1
Construction	
Material	Stainless steel
Weight	11 110 g
Nominal value	1 g – 5 kg
Quantity	16 pcs

## Compatible with (Additional Fee)

Tweezers

Brushes