



More information on the website  
[mirror.radwag.com/de/info,w1,X6X](http://mirror.radwag.com/de/info,w1,X6X)

# Prüfgewicht der E2 Klasse - Knopfgewicht - 20 g

OP-103-0042



The drawings, photos and graphics used are for illustrative purposes only.

## Technische Daten

Messtechnische Parameter	
OIML-Klasse	E2
Konstruktion	
Material	Edelstahl
Nominalwert	20 g

## Kompatibel mit (Additional Fee)

Pinzetten

Pinself