



More information on the website  
[mirror.radwag.com/en/info,w1,XE2](https://mirror.radwag.com/en/info,w1,XE2)

# F1 Mass Standard - Knob Weights With Adjustment Chamber, Set (1 mg - 200 g), Wooden Box

OK-201-0055



The drawings, photos and graphics used are for illustrative purposes only.

## Datasheet

Metrological parameters	
OIML Class	F1
Construction	
Material	Stainless steel
Weight	611.11 g
Nominal value	1 mg – 200 g
Quantity	23 pcs

## Compatible with (Additional Fee)

Tweezers

Brushes