



More information on the website
mirror.radwag.com/en/info,w1,XF7

Knob weights, Set (1 mg - 5 kg), Wooden Box

OK-101-0045



The drawings, photos and graphics used are for illustrative purposes only.

Datasheet

| Metrological parameters | |
|-------------------------|-----------------|
| OIML Class | E2 |
| Construction | |
| Material | Stainless steel |
| Weight | 11 111.11 g |
| Nominal value | 1 mg – 5 kg |
| Quantity | 28 pcs |

Compatible with (Additional Fee)

Tweezers

Brushes