



More information on the website
mirror.radwag.com/en/info,w1,XM1

PT0443 Cable



The drawings, photos and graphics used are for illustrative purposes only.

Overview

PT0443 Cable is used to connect the 5Y series balance or the PUE CY10 terminal to a computer. Its standard length is 1,8 m.

Optionally, the cable is used to connect the 5Y series balance or the PUE CY10 terminal to:

AS X2 and AS R series balances PS X2 i PS R series balances PM C32 series balances WLC and WTC series balances, and C315 series scales MAS and MPS series weighing modules

Compatible with

H315 HR Waterproof Platform Scale
 MAS1 Modules
 XA 5Y.M.A.P Microbalance for pipettes calibration
 MYA 5Y.P Microbalance for pipettes calibration
 Tunnel Detectors
 C315 Platform scale
 5Y.PM Precision Balance
 MYA 5Y.F Microbalance for Filters
 AS X7 Analytical Balance

PS X7 Precision Balance
 XA 5Y.A Analytical Balance
 C315.4P pallet scales
 UYA 5Y Ultra-Microbalance
 AS 5Y Analytical Balance
 PS 5Y Balances for PGC
 PS R2.H Precision Balance
 WTC Precision Balance
 CY10 Precision Balance

PS X2 Precision Balance
AS X2 PLUS Analytical Balances
PS X7 Balances for PGC
C315.D Baby Scale
XA 5Y.M.A Microbalance
WLC Precision Balance
PS R2 Precision Balance
C315.4P2 Beam Scales
H315 4H Stainless Steel Platform Scale
H315 4P2.H Stainless steel Beam Scales
WLC X7 Precision Balance
CY10 Scales for PGC
PS X2 Balances for PGC
H315 H Waterproof Platform Scale
PUE CY10 Weighing Terminal
CY10 Multifunctional Scale
MAS Modules
H315.4 H/Z Stainless Steel Platform Scale, pit version
MPS Modules
H315 4N.H Stainless Steel Ramp Scale
C315.4 Platform Scale
C32.PM Precision Balance
C315.O Personal Scale
H315 4N.H.LD Stainless Steel Ramp Scale

H315 4C Platform Scale
C315.OW Personal Scale
XA 5Y.F Analytical Balance for Filters
H315 2K Overhead Track Scales
H315 4P.H Stainless Steel Pallet Scales
WLC C/2 Precision Balance
MYA 5Y Microbalance
H315 4I.S Livestock Scales
PS R1 Precision Balance
XA 5Y Analytical Balance
PS 5Y Precision Balance
Postal Scales for Packages
PMV 5Y Moisture Analyzer
H315 4I.H Stainless Steel Livestock Scales
XA 5Y.M Microbalance
PUE C315 Weighing Terminal
AS R1 PLUS Analytical Balance
UYA 5Y.F Ultra-Microbalance for Filters
C315.4B Bed Scale
PUE H315 Weighing Terminal
C315.K Chair Scale
XA 5Y.M.A.S Microbalance for stents weighing
H315 Stainless Steel Overhead Track Scales
C315.8B Ramps Bed Scale