



More information on the website  
[mirror.radwag.com/en/info,w1,XM1](http://mirror.radwag.com/en/info,w1,XM1)

## PT0443 Cable

WX-013-0229



The drawings, photos and graphics used are for illustrative purposes only.

## Overview

PT0443 Cable is used to connect the 5Y series balance or the PUE CY10 terminal to a computer. Its standard length is 1,8 m.

Optionally, the cable is used to connect the 5Y series balance or the PUE CY10 terminal to:

AS X2 and AS R series balances PS X2 i PS R series balances PM C32 series balances WLC and WTC series balances, and C315 series scales MAS and MPS series weighing modules

## Compatible with (Additional Fee)

H315 HR Waterproof Platform Scale  
MAS1 Modules  
XA 5Y.M.A.P Microbalance for pipettes calibration  
MYA 5Y.P Microbalance for pipettes calibration  
Tunnel Detectors  
C315 Platform scale  
5Y.PM Precision Balance  
MYA 5Y.F Microbalance for Filters  
AS X7 Analytical Balance

PS X7 Precision Balance  
XA 5Y.A Analytical Balance  
C315.4P pallet scales  
UYA 5Y Ultra-Microbalance  
AS 5Y Analytical Balance  
PS 5Y Balances for PGC  
PS R2.H Precision Balance  
WTC Precision Balance  
CY10 Precision Balance

PS X2 Precision Balance  
AS X2 PLUS Analytical Balances  
PS X7 Balances for PGC  
C315.D Baby Scale  
XA 5Y.M.A Microbalance  
WLC Precision Balance  
PS R2 Precision Balance  
C315.4P2 Beam Scales  
H315 4H Stainless Steel Platform Scale  
H315 4P2.H Stainless steel Beam Scales  
WLC X7 Precision Balance  
CY10 Scales for PGC  
PS X2 Balances for PGC  
H315 H Waterproof Platform Scale  
PUE CY10 Weighing Terminal  
CY10 Multifunctional Scale  
MAS Modules  
H315.4 H/Z Stainless Steel Platform Scale, pit version  
MPS Modules  
H315 4N.H Stainless Steel Ramp Scale  
C315.4 Platform Scale  
C32.PM Precision Balance  
C315.O Personal Scale  
H315 4N.H.LD Stainless Steel Ramp Scale

H315 4C Platform Scale  
C315.OW Personal Scale  
XA 5Y.F Analytical Balance for Filters  
H315 2K Overhead Track Scales  
H315 4P.H Stainless Steel Pallet Scales  
WLC C/2 Precision Balance  
MYA 5Y Microbalance  
H315 4I.S Livestock Scales  
PS R1 Precision Balance  
XA 5Y Analytical Balance  
PS 5Y Precision Balance  
Postal Scales for Packages  
PMV 5Y Moisture Analyzer  
H315 4I.H Stainless Steel Livestock Scales  
XA 5Y.M Microbalance  
PUE C315 Weighing Terminal  
AS R1 PLUS Analytical Balance  
UYA 5Y.F Ultra-Microbalance for Filters  
C315.4B Bed Scale  
PUE H315 Weighing Terminal  
C315.K Chair Scale  
XA 5Y.M.A.S Microbalance for stents weighing  
H315 Stainless Steel Overhead Track Scales  
C315.8B Ramps Bed Scale