



# PT0443 Cable

WX-013-0229

More information on the website  
[mirror.radwag.com/us/info,w1,XM1](http://mirror.radwag.com/us/info,w1,XM1)



The drawings, photos and graphics used are for illustrative purposes only.

## Overview

PT0443 Cable is used to connect the 5Y series balance or the PUE CY10 terminal to a computer. Its standard length is 1,8 m.

Optionally, the cable is used to connect the 5Y series balance or the PUE CY10 terminal to:

AS X2 and AS R series balances PS X2 i PS R series balances PM C32 series balances WLC and WTC series balances, and C315 series scales MAS and MPS series weighing modules

## Compatible with (Additional Fee)

H315 HR Waterproof Platform Scale  
 XA 5Y.M.A.P Microbalance for pipettes calibration  
 MYA 5Y.P Microbalance for pipettes calibration  
 Tunnel detectors  
 C315 Platform scale  
 5Y.PM Precision balance  
 MYA 5Y.F Microbalance for Filters  
 AS X7 Analytical Balance  
 PS X7 Balances for PGC  
 C315.D Baby scale

AS 5Y Analytical Balance  
 PS 5Y Precision Balance  
 PS 5Y Balances for PGC  
 PS R2.H Precision Balance  
 WTC Precision Balance  
 CY10 Precision balance  
 H315 4C Platform Scale  
 C315.OW Personal scale  
 Analytical Balance for Filters XA 5Y.F  
 H315 2K Overhead Track Scales

XA 5Y.M.A Microbalance  
WLC Precision Balance  
PS R2 Precision Balance  
C315.4P2 Beam Scales  
H315 4H Stainless Steel Platform Scale  
H315 4P2.H Stainless steel Beam Scales  
WLC X7 Precision balance  
CY10 Scales for PGC  
H315 H Waterproof Platform Scale  
PUE CY10 Weighing Terminal  
CY10 Multifunctional Scale  
H315.4 H/Z Stainless Steel Platform Scales, pit version  
H315 4N.H Stainless Steel Ramp Scale  
C315.4 Platform Scale  
C32.PM Precision Balance  
C315.0 Personal scale  
H315 4N.H.LD Stainless Steel Ramp Scale  
XA 5Y.A Analytical Balance  
C315.4P pallet scales  
UYA 5Y Ultra-Microbalance

PS X7 Precision Balance  
H315 4P.H Stainless Steel Pallet Scales  
WLC C/2 Precision balance  
MYA 5Y Microbalance  
PS R1 Precision Balance  
XA 5Y Analytical Balance  
Postal scales for packages  
PMV 5Y Moisture Analyzer  
XA 5Y.M Microbalance  
PUE C315 Terminal  
AS R1 PLUS Analytical Balance  
UYA 5Y.F Ultra-Microbalance for Filters  
C315.4B Bed scale  
PUE H315 Terminal  
C315.K Chair scale  
XA 5Y.M.A.S Microbalance for stents weighing  
MAS.1.Y Modules  
H315 2K Stainless steel overhead track scales  
C315.8B Ramps bed scale