



More information on the website
mirror.radwag.com/us/info,w1,XM1

PT0443 Cable

WX-013-0229



The drawings, photos and graphics used are for illustrative purposes only.

Overview

PT0443 Cable is used to connect the 5Y series balance or the PUE CY10 terminal to a computer. Its standard length is 1,8 m.

Optionally, the cable is used to connect the 5Y series balance or the PUE CY10 terminal to:

AS X2 and AS R series balances PS X2 i PS R series balances PM C32 series balances WLC and WTC series balances, and C315 series scales MAS and MPS series weighing modules

Compatible with (Additional Fee)

H315 HR Waterproof Platform Scale
XA 5Y.M.A.P Microbalance for pipettes calibration
MYA 5Y.P Microbalance for pipettes calibration
Tunnel detectors
C315 Platform scale
5Y.PM Precision balance
MYA 5Y.F Microbalance for Filters
AS X7 Analytical Balance
PS X7 Balances for PGC
C315.D Baby scale

UYA 5Y Ultra-Microbalance
AS 5Y Analytical Balance
PS 5Y Balances for PGC
PS R2.H Precision Balance
WTC Precision Balance
CY10 Precision balance
H315 4C Platform Scale
C315.OW Personal scale
Analytical Balance for Filters XA 5Y.F
H315 2K Overhead Track Scales

XA 5Y.M.A Microbalance
WLC Precision Balance
PS R2 Precision Balance
C315.4P2 Beam Scales
H315 4H Stainless Steel Platform Scale
H315 4P2.H Stainless steel Beam Scales
WLC X7 Precision balance
CY10 Scales for PGC
H315 H Waterproof Platform Scale
PUE CY10 Weighing Terminal
CY10 Multifunctional Scale
H315.4 H/Z Stainless Steel Platform Scales, pit version
H315 4N.H Stainless Steel Ramp Scale
C315.4 Platform Scale
C32.PM Precision Balance
C315.0 Personal scale
H315 4N.H.LD Stainless Steel Ramp Scale
PS X7 Precision Balance
XA 5Y.A Analytical Balance
C315.4P pallet scales

H315 4P.H Stainless Steel Pallet Scales
WLC C/2 Precision balance
MYA 5Y Microbalance
PS R1 Precision Balance
XA 5Y Analytical Balance
PS 5Y Precision Balance
Postal scales for packages
MAS1 Modules
PMV 5Y Moisture Analyzer
XA 5Y.M Microbalance
PUE C315 Terminal
AS R1 PLUS Analytical Balance
UYA 5Y.F Ultra-Microbalance for Filters
C315.4B Bed scale
PUE H315 Terminal
C315.K Chair scale
XA 5Y.M.A.S Microbalance for stents weighing
H315 Stainless steel overhead track scales
C315.8B Ramps bed scale