



More information on the website
mirror.radwag.com/en/info,w1,XRY

KR-01 Converter

WX-004-0021



The drawings, photos and graphics used are for illustrative purposes only.

Datasheet

Physical parameters

Packaging dimensions	220×190×90 mm
Net weight	0.4 kg
Gross weight	0.8 kg

Compatible with (Additional Fee)

MA X7.IC.A Moisture Analyzer
 MA 5Y Moisture Analyzer
 Tunnel Detectors
 C315 Platform scale
 MA X2.IC.A Moisture Analyzers
 AS X7 Analytical Balance
 PS X2 Precision Balance
 AS X2 PLUS Analytical Balances
 PS X7 Balances for PGC
 C315.D Baby Scale
 WLC Precision Balance
 C315.4P2 Beam Scales
 WLC X7 Precision Balance
 PS X2 Balances for PGC

C315.0 Personal Scale
 PS X7 Precision Balance
 C315.4P pallet scales
 MA 5Y.IC Moisture Analyzer
 AS 5Y Analytical Balance
 PS 5Y Balances for PGC
 H315 4C Platform Scale
 C315.0W Personal Scale
 WLC C/2 Precision Balance
 PS 5Y Precision Balance
 Postal Scales for Packages
 PUE C315 Weighing Terminal
 C315.4B Bed Scale
 PUE H315 Weighing Terminal

MA X2.A Moisture Analyzers
C315.4 Platform Scale
MA X7.A Moisture Analyzer

C315.K Chair Scale
C315.8B Ramps Bed Scale

Device dimensions

