



More information on the website
mirror.radwag.com/us/info,w1,XRY

KR-01 Converter

WX-004-0021



The drawings, photos and graphics used are for illustrative purposes only.

Datasheet

Physical parameters

Packaging dimensions W x D x H	220×190×90 mm
Net weight	0,4 kg
Gross weight	0,8 kg

Compatible with (Additional Fee)

MA X7.IC.A Moisture Analyzer
 MA 5Y Moisture Analyzer
 Tunnel detectors
 C315 Platform scale
 AS X7 Analytical Balance
 PS X7 Balances for PGC
 C315.D Baby scale
 WLC Precision Balance
 C315.4P2 Beam Scales
 WLC X7 Precision balance
 C315.4 Platform Scale
 MA X7.A Moisture Analyzer
 C315.O Personal scale
 C315.4P pallet scales

MA 5Y.IC Moisture Analyzer
 AS 5Y Analytical Balance
 PS 5Y Precision Balance
 PS 5Y Balances for PGC
 H315 4C Platform Scale
 C315.OW Personal scale
 PS X7 Precision Balance
 WLC C/2 Precision balance
 Postal scales for packages
 PUE C315 Terminal
 C315.4B Bed scale
 PUE H315 Terminal
 C315.K Chair scale
 C315.8B Ramps bed scale

Device dimensions W x D x H

