



More information on the website  
[mirror.radwag.com/en/info,w1,XYU](http://mirror.radwag.com/en/info,w1,XYU)

# Knob weights, Set (1 g - 5 kg), Plastic Box

OK-102-0029



The drawings, photos and graphics used are for illustrative purposes only.

## Datasheet

Metrological parameters	
OIML Class	E2
Construction	
Material	Stainless steel
Weight	11 110 g
Nominal value	1 g – 5 kg
Quantity	16 pcs

## Compatible with (Additional Fee)

Tweezers

Brushes