



More information on the website  
[mirror.radwag.com/en/info,w1,Y5W](https://mirror.radwag.com/en/info,w1,Y5W)

# 105 mm carbon fiber tweezers for mass standards up to 50 g

OA-105-0001



The drawings, photos and graphics used are for illustrative purposes only.

## Datasheet

Construction	
Material	carbon fiber
Physical parameters	
Length	105 mm

## Compatible with (Additional Fee)

E1 Mass Standard - knob weights, sets (1 mg - 10 kg)

F1 Mass Standard - stainless steel wires

F1 Mass Standard - knob weights with adjustment chamber, sets

F2 Mass Standard - stainless steel wires

F2 Mass Standard - Carrier weights

E2 Mass Standard - stainless steel wires

E2 Mass Standard -

F1 Mass Standard -

knob weights with adjustment chamber

E2 Mass Standard - flat polygonal sheets

E2 Mass Standard -

knob weights

F1 Mass Standard - flat polygonal sheets, sets

F1 Mass Standard -

knob weights without adjustment chamber

E1 Mass Standard -

stainless steel wires

F2 Mass Standard -

knob weights with adjustment chamber

F2 Mass Standard - flat polygonal sheets, sets

knob weights, sets (1 mg - 10 kg)  
F2 Mass Standard -  
knob weights with adjustment chamber, sets  
E2 Mass Standard -  
cylindrical weights, sets (1 g - 10 kg)  
F2 Mass Standard -  
Slotted weights  
M1 Mass Standard - flat polygonal sheets  
F2 Mass Standard - flat polygonal sheets  
F1 Mass Standard -  
nob weights without adjustment chamber, set  
F1 Mass Standard - flat polygonal sheets  
E1 Mass Standard -  
stainless steel wires, sets (1 mg - 500 mg)  
F1 Mass Standard -  
cylindrical weights  
E2 Mass Standard -  
stainless steel wires, sets (1 mg - 500 mg)

M1 Mass Standard -  
knob weights  
E1 Mass Standard - flat polygonal sheets  
E1 Mass Standard -  
knob weights  
E2 Mass Standard - flat polygonal sheets, sets (1 mg - 500 mg)  
E1 Mass Standard - flat polygonal sheets, sets (1 mg - 500 mg)  
F2 Mass Standard -  
Stainless steel wires, sets  
F1 Mass Standard -  
stainless steel wires, sets  
M1 Mass Standard - flat polygonal sheets, sets